

Muscle Activity in Ipsi-and Contralateral limbs in Total Knee Arthroplasty Patient during Level and Sloped Walking

Mario Lamontagne, PhD^{1,2,3} | Erik Kowalski, MSc¹ | Clara Philips² | Geoffrey Dervin, MD³ mlamon@uottawa.ca | e.kowalski@uottawa.ca | cphili081@uottawa.ca | gdervin@toh.ca

Background

- Total knee arthroplasty (TKA) patients with medial pivot (MP) implants provide more similar muscle activity patterns to healthy knees than posterior stabilized (PS) prostheses.^{1,2}
- However, it is unknown how the prosthesis design would affect the muscle activity of the nonoperated limb.

Tables

Table 1: Participant demographics

Group	Number	Age (years)	BMI (kg/m²)	Post-Op (months)
MP	8	59.9 ± 5.5	29.0 ± 4.5	10.8 ± 3.0
PS	6	68.3 ± 3.6	29.8 ± 3.0	9.8 ± 1.7
СТ	9	65.7 ± 5.4	27.1 ± 5.0	N/A

CIHR IRSC

Objective

• The purpose of this study was to compare lower limb muscle activity in patients who underwent a total knee arthroplasty (TKA) with a medial pivot (MP) implant to healthy controls (CTRL) during a

Results

• PS group had significantly higher iEMG in the quadriceps and hamstrings of the non-operated limb during the level walking task, whereas MP group had no significant differences for the same muscles

stair ascent task.

Methods

- 14 participants required a unilateral total knee arthroplasty and were randomly assigned to either a MP (MicroPort Orthopaedics Inc.) or a PS (Zimmer Biomet) implant and were compared to 9 healthy controls (CT) (Table 1)
- Muscle activities were collectured using wireless electromyography (EMG) probes on the semimembranosus and biceps femoris (hamstings); vastus medialis and vastus lateralis (quadriceps); medial and lateral heads of of the gastrocnemius (gastrocnemii) of both limbs
- EMG data was recorded as participants completed level and sloped (±9°) walking
- Total muscle activity (iEMG) was extrapolated and normalized to a maximal voluntary isometric contraction
- Significant differences between operated and non-operated limbs in TKA group and left and right limbs in the CTRL group was obtained using a paired t-test ($\alpha < 0.05$)







• MP group had significantly higher iEMG in the hamstrings and gastrocnemii of the operated limb than the non-operated limb during the inclined walking task

Discussion

- Previous studies have found that TKA patients have higher and longer muscle activity in the nonoperated limb compared to the operated limb during walking^{3,4}
- PS group resulted in greater iEMG on non-operated limb compared to operated limb during walking as was previously found¹
 - This may indicate a compensatory gait originating from muscle adaptations resulting of many years of waiting for the knee replacement
- MP group had better inter-limb muscle symmetry for level and declined walking, which may indicate a greater ability to load both limbs equally, reducing overload on the contralateral limb
- Differences in implant designs between the MP and PS implant may contribute to differences in muscle activation
- Reduced muscle activation on the non-operated limb may reduce joint loading, consequently limiting joint degeneration
- Conformity in implant designs showed less kinematic and kinetic variability⁵, consequently lowering joint loading
- Therefore, post-surgery rehabilitation should aim to strengthen muscles in the operated limb so



Figure : Total muscle activity for the MP, PS and CT groups during the three walking conditions. * represents significant differences between operated and non-operated limbs in the MP and PS group, or between left and right limbs of the CT group







Acknowledgements

The authors wish to acknowledge the funding contributions from the Canadian Institutes of Health Research, the Natural Sciences and Engineering Research Council of Canada and MicroPort Orthopedics Inc.



Schmidt et al (2003) *Clinical Orthopaedics and Related Research* Omori et al (2009) *Journal of Orthopedic Science* Bendetti et al (2003) *Clinical Biomechanics*

Stevens-Lapsley et al (2010) *Clinical Orthopaedics and Related Research* Ardestani et al (2015) *Medical Engineering & Physics*

COA-ACO

2017 OTTAWA

JUNE 15-18 JUIN

